

Specifications of L-0189

●MICROSCOPE

Type:	Galilean-Type Binocular Microscope
Magnification Changer:	5 Steps by drum rotation
Eyepieces:	12.5X
Total Magnifications:	6, 10, 16, 25, 40X
Field of View:	38,24,15,9,6mm
Pupillary Adjustment:	48-78mm
Focus Distance:	100 mm
Diopter Adjustment:	+5D- -6D

●SLIT ILLUMINATION

Slit Width:	0-10mm
Aperture Diameter:	0.2, 1, 3, 4, 6, 10mm
Slit Length:	0-10mm
Slit Rotation:	±90°
Filters:	Cobalt blue, Heat absorption, Red-free, 1/2ND
Slit Inclination:	5° 10° 15° 20°
Lamp:	12V/50W Halogen Lamp(L0160H1)

●BASE

Vertical Movement:	30mm
Longitudinal Movement:	70mm
Lateral Movement:	115mm
Fine Base Movement:	±10mm

●CHINREST

Vertical Movement:	75mm
Fixation Target:	Tungsten lamp 6V 0.15A
Input Voltage:	AC. 100, 110, 220, 240V 50/60Hz
Power Consumption:	50VA

●DIMENSION & WEIGHT

Weight:	20.5kgs
Dimension:	550(W)x350(D)x500(H)



Certificate:
09 105 89060

Inami & Co., Ltd.

L-0189

5 STEPS MAGNIFICATION SLIT LAMP

5 STEPS MAGNIFICATION SLIT LAMP L-0189



L-0240
ZOOM MAGNIFICATION SLIT LAMP



Certificate:
09 105 89060



●Specifications subject to change without prior notice.



Inami & Co., Ltd.

BIJYOGI RESEARCH INSTITUTE
No. 3-14-13, Bijyogi, Toda-City,
Saitama-Pref, 335-0031, Japan
Tel. 81-48-422-0173 Fax. 81-48-421-3838
Website : <http://www.inami.co.jp>
E-Mail : ibri@inami.co.jp
®TRADE MARK OF INAMI & CO., LTD.

HEAD OFFICE
No. 24-2, Hongo 3-Chome,
Bunkyo-ku, Tokyo 113-0033, Japan
Tel. 81-3-3814-5916 Fax. 81-3-5684-2126
Website : <http://www.inami.co.jp>
E-Mail : corpinfo@inami.co.jp
®TRADE MARK OF INAMI & CO., LTD.





MS-4
5 STEPS MAGNIFICATION SYSTEM
(6, 10, 16, 25, 40X, OBJECTIVE LENS
f=100mm)



OCL-004
EYE PIECES 12.5X



L-0542-C
ASSISTANT BINOCULAR MICROSCOPE



L-0541
BEAM SPLITTER



L-0544
TV CAMERA ADAPTER(f=100mm)

The L-0189 5 STEPS MAGNIFICATION SLIT LAMP provides practitioners with the utmost satisfaction thanks to its excellent image clearness and fine design.



YF-57
SPECIAL ADAPTOR BUILT-IN
YELLOW FILTER



VARIOUS KINDS OF BUILT-IN FILTERS:
COBALT BLUE, RED-FREE, 50%ND AND
HEAT ABSORPTION FILTERS

